# *Mini LED*

Auto.Focus

---

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# ENGLISH

# SUMMARY

1	SAFETY INSTRUCTIONS	6	10	TECHNICAL SERVICE	10
2	PRESENTATION	7	11	SPECIFICATIONS	11
3	INSTALLATION	7	12	ELECTROMAGNETIC COMPATIBILITY	11
4	DESCRIPTION	7	13	DISPOSAL AND RECYCLING	14
5	MENUS	8	14	LIABILITY	14
6	ROUTINE USE	8	15	REGULATIONS	14
7	CONFIGURATION	9	16	SYMBOLS AND ABBREVIATIONS	14
8	MAINTENANCE	9	17	CUSTOMER RELATIONS	74
9	TROUBLESHOOTING	10			

# 1 - SAFETY INSTRUCTIONS

**Warning:** Before using your **Mini L.E.D. AutoFocus**, please read the whole of this Safety Instructions section carefully. The manufacturer's warranty is applicable only if these instructions concerning the operation and the safety of the unit have been properly applied. All these safety measures assume a good level of knowledge of dental techniques and photo-polymerization, and the more specific knowledge regarding the use of the **Mini L.E.D. AutoFocus** contained in this Operating Manual.

- 1) This unit must be used in strict accordance with the instructions contained in this Operating Manual. The manufacturer declines any liability if these rules are not followed or if the unit is used for any other application.
- 2) Before connecting the unit, make sure that the power voltage is compatible with the voltage marked on the power adapter used to charge the curing light. A different voltage would damage the unit and could result in injury to the patient and/or the operator.
- 3) The light radiation produced by this type of unit can be dangerous and must not be pointed at the eyes, even if the practitioner or the patient is wearing suitable protective eyewear (class 2M laser). The light produced by this equipment must be directed only at the zone to be treated in the oral cavity.
- 4) Any condensation inside electrical equipment can be a hazard. If the curing light has to be moved from a cool place to a warm place, do not use it immediately: let it reach the new room temperature first.
- 5) Do not replace the battery during use.
- 6) Comply with all the instructions given in this Manual. Do not insert or attempt to insert metal objects into the unit or its base station, as this could cause an electrical discharge.
- 7) The manufacturer declines all liability if damaged parts or accessories are not replaced exclusively by the manufacturer's products. In particular, the use of other light guides, power adapters or batteries could be dangerous for the patient and the operator.
- 8) In the case of a fault, unplug the base station immediately and make sure that no one can use the curing light before it is checked by the manufacturer or dealer. The fault may be the result of non-compliance with the safety rules or of a technical problem in the unit.
- 9) Do not use the unit near a heat source. Do not use solvents, detergents or flammable products, which may damage the unit or cause short-circuits.
- 10) Only the unit manufacturer or dealer is qualified to carry out repairs.
- 11) All lights used for photo-polymerization, including the **Mini L.E.D. AutoFocus**, must not be used on persons who are suffering or have suffered from photo-biological reactions (including urticaria solaris or erythropoietic protoporphyrina) or who are receiving treatment including photosensitizing medications (including methoxsalens or chlorotetracycline).
- 12) Persons (practitioners or patients) who have previously suffered from a retina or crystalline lens condition or who have undergone eye surgery, in particular for cataracts, must see their ophthalmologist before using the **Mini L.E.D. AutoFocus**. Even if the eye specialist agrees, prudence is strongly recommended, because the intensity of the light could cause accidents. Class 2M protective eyewear, suitable for use with units emitting radiation at wavelengths between 420 and 480 nm, should be worn at all times.
- 13) **Mini L.E.D. AutoFocus** must not be used if the patient and/or the operator has a cardiac stimulator or any other active implant (e.g. a cochlear implant).
- 14) The device is not designed to withstand shocks delivered by an electric defibrillator.
- 15) Any variation of the power supply voltage or the magnetic field outside the limits in force could switch the curing light to automatic operation or interfere with its operation.
- 16) Electromagnetic interference: the unit complies with the applicable standards (IEC 60 601-1-2) for emission (radiated electrical field and interference with the supply voltage) and for immunity (protection against electrical fields, fast supply transients in bursts, electrostatic discharge and supply surges).
- 17) Your dealer carries stocks of all spare parts that you might need. Used batteries must be returned to the dealer.
- 18) For use only by dental professionals.
- 19) When you need to transport the unit, unscrew the battery and protect the light guide against impacts.
- 20) For long-term storage we recommend that you remove the battery from the unit body so that it does not discharge slowly and risk damage.

- 21) When your unit reaches the end of its life, we recommend that you contact your dental equipment dealer or, failing this, one of the ACTEON sites (listed at the end of this manual) to find out what procedure is advised.

## 2 - PRESENTATION

The Mini L.E.D. AutoFocus emits visible blue light at wavelengths between 420 and 480nm for photopolymerization of dental materials. After checking the packaging (please keep it) for damage, open the box, which should contain the following items:

- the handpiece, fitted with a protective cap, incorporating the light-emitting diode and the back-lit LCD screen (3),
- the base station for recharging the unit, with a built-in radiometer,
- the power adapter,
- the lithium ion battery (7),
- the sterilizable light guide (1),  
the flexible filter shield,
- the Operating Manual and other documents.

## 3 - INSTALLATION

- 1) Before use, sterilize the light guide and disinfect the base station and the unit (refer to the Maintenance section).
- 2) Screw the battery onto the handpiece and insert the sterilized light guide into the unit. Make sure that the light guide is inserted correctly (confirmation by an audible click).
- 3) Place the base station on a stable surface. Check that the voltage marked on the power adapter is the same as your power supply voltage, then plug in the adapter. The translucent holder flashes three times (red, green, blue) and emits a beep at power-on.
- 4) Place the curing light on its rest. Two beeps confirm that it is positioned correctly. A complete battery charging cycle lasts 3 hours.
- 5) The translucent holder lights up in blue and starts flashing, indicating that the battery is charging. When the battery is fully charged, the holder stops flashing and the blue light stays on.

## 4 - DESCRIPTION

Now that the Mini L.E.D. AutoFocus is fully charged,

you need to familiarize yourself with its controls and how to use it.  
(See Figure 1 and 2, page 76-77)

### Handpiece:

#### Buttons

The handpiece has three buttons:

- an ON/OFF button (2) to start or stop the polymerization cycle


ON/OFF BUTTON


- selector buttons (5) & (6) for selecting the type of menu and the unit configuration

Configuration button

#### LCD display (3):

The LCD display shows, from top to bottom:

- readout of light emission time for polymerization (in seconds),
- three pictograms indicating the selected polymerization mode,
- symbols AF and LA, which appear if the AutoFocus and Laser Target are activated,
- diameter of the selected light guide: 7.5mm (larger circle) or 5.5mm.

The table below tells you which polymerization mode is indicated by each pictogram.

Pictogram	Mode
	Fast curing
	Pulse curing
	Stepped curing

#### Indicator light:

The handpiece has a two-color indicator light (4) located in the middle of the double selector button (5 &

6). The table below gives the operating status of the curing light indicated by the indicator light:

Light color	Status
Green	Normal mode
Red	Battery low
Red flashing	Thermal protection

#### **Base station:**

The base station is used to recharge the battery. The translucent handpiece rest (8) also acts as an indicator. The station incorporates a radiometer, which you can use to check whether your curing light fitted with its light guide is powerful enough to ensure satisfactory polymerization. The radiometer window is located on the front of the rest (9).

(See Figure 2 of the page 77)

#### **Radiometer:**

To use the radiometer, proceed as follows:

- set the curing light to Fast curing mode,
- check that the standard light guide is intact and bears no traces of composite material, then insert the guide into the handpiece,
- place the end of the light guide flat against the radiometer window (9) and switch on your **Mini LED AutoFocus**,
- note the color with which the handpiece holder (8) lights up, depending on the measured power:

Color	Status
Green	Power ok
Red	Refer to section 9: Troubleshooting

When the handpiece is placed on its rest for charging, the translucent rest (8) lights up as follows:

Color	Status
Blue flashing	Battery charging
Blue steady	Battery charged

## **5 - MENUS**

The **Mini LED AutoFocus** has three different menus. Press M on the double selector button (5) to select the menu you want to use.

"*Fast curing*" mode: When the handpiece is removed

from the base unit to be used for the first time, the curing light is in Fast curing mode. The pictogram on the LCD display indicates the mode. In Fast curing mode the curing light power is:

- 2,200 mW/cm<sup>2</sup> (± 10%) with the standard 7.5mm light guide,
- 3,000 mW/cm<sup>2</sup> (± 10%) with the optional Ø 5.5 mm light guide.

"*Pulse curing*" mode: Press M on the double selector button (5) to select Pulse curing. This menu gives the same power as the curing light delivers in Fast curing mode, but in 10 consecutive flashes with a 250 ms rest period between flashes.

"*Stepped curing*" mode: Press M on the double selector button (5) to select Stepped curing. In Stepped curing mode (polymerization similar to the Soft menu with a halogen curing light) you have a ramp-up to full power over 10 seconds then full power for 10 seconds, at the level determined by the selected light guide:

- from 0 to 2,200mW/cm<sup>2</sup> with the standard 7.5mm light guide,
- from 0 to 3,000mW/cm<sup>2</sup> with the optional 5.5mm light guide.

## **6 - ROUTINE USE**

The **Mini L.E.D. AutoFocus** has two particularly innovative features enabling you to optimize the polymerization of your composites:

- a Laser Target function to aim the most powerful part of your beam at a particular position on your photo-polymerizable material,
- an AutoFocus function that automatically adjusts the light emission time to optimize the polymerization energy delivered to your material. AutoFocus enables you to obtain constant polymerization quality to a depth of 2mm, correcting the positioning variations there may be between the light guide tip and the clinical site.

1) The **Mini L.E.D. AutoFocus** must normally be placed on its rest. The first time the handpiece is screwed onto the battery, the Fast curing mode is displayed. For more details refer to section 5.


2) Select the menu to use by pressing button M (5), then select the light guide used by pressing the configuration button (6). The selected light guide is symbolized by a filled circle displayed at the bottom of the LCD screen. The diameter shown on the display

must correspond to that marked on the light guide fitting. This is essential for proper operation of the AutoFocus function. The larger filled circle pictogram on the left indicates the standard Ø7.5mm light guide, while the smaller filled circle pictogram on the right indicates the optional Ø5.5mm light guide.

- 3) Bring the light guide close to the material surface, without touching it as that would risk reducing polymerization quality. The use of a transparent strip is advised.
- 4) Press the ON/OFF button to start the polymerization cycle. This action is confirmed by an audible beep.
- 5) Press and hold the ON/OFF button to produce a red aiming circle. When you release the button, the aiming circle disappears and the curing light is switched on. This function allows you to view the zone where the power is maximal.

- 6) The AutoFocus function needs a one-second adjustment time to determine the polymerization parameters. The preprogrammed polymerization time of the selected menu varies according to the position of the light guide with respect to the material surface. The calculated time in seconds is displayed on the LCD screen.

- 7) When the adjustment time has elapsed, the cycle is started automatically. A countdown on the LCD screen shows you how much time remains, and every five seconds a beep informs you of the progress of the current cycle. In AutoFocus mode, when the polymerization cycle is completed, the time that was calculated and applied is displayed.
- 8) You can press the ON/OFF button at any time to interrupt the current polymerization cycle.
- 9) At the end of the cycle, another beep informs you that the cycle has been completed successfully.
- 10) Note that if no button has been pressed for 30 seconds the LCD screen backlighting switches off.

11) After three minutes without use, the unit switches to standby mode (low consumption). Press any button to reactivate the light.

## 7 - CONFIGURATION

The AutoFocus and Laser Target function settings can be defined by the user: you can decide whether to activate them or not.

### AutoFocus configuration

The AutoFocus function is active only in Fast curing and Pulse curing modes. The range of measurement of the AutoFocus function is 9mm.

To change the status of the AF mode, press and hold the configuration button (6) of the double button for 2.5 seconds.

The screen below is displayed:

When the function is active, the symbol AF is displayed below the menu pictograms.

When the Laser Target and AutoFocus functions are activated at the same time, the Target ring appears then disappears when the ON/OFF button is released. In this case, releasing the button starts the AutoFocus function measurement.

### Laser Target configuration

By default, the Mini L.E.D. AutoFocus is configured with Laser Target active.

To activate or deactivate Laser Target, press and hold the menu button (5) for 2.5 seconds.

The screen below is displayed:

If the pictogram LA is visible, Target is active.

## 8 - MAINTENANCE

Warning: For everyday maintenance of the Mini L.E.D. AutoFocus, make sure that the unit is not on its rest.

Disconnect the rest from the Mini L.E.D. AutoFocus before using a disinfectant.

1) The Mini L.E.D. AutoFocus and its accessories can be cleaned and disinfected with ready-to-use cleaning and disinfecting wipes impregnated with alcohol, an amphoteric and a biguanide (refer to manufacturer's instructions) for at least two minutes. Leave the product to take effect for at least 15 minutes. Use wipes that have CE marking or that are in compliance with any standard required by the national regulations.

Wipe dry using a disposable clean non-woven cloth until there are no traces of liquid.

Packaging: Keep sterilizable accessories in single-use sterilization sleeves or bags in compliance with the specifications defined in standard EN ISO 11607-1 or in any equivalent standard required by the national regulations.

Sterilize the light guide and the flexible optical shield individually in an autoclave before each patient, according to the following parameters:

- autoclave, compliant with standard EN 13060, class B
- sterilization temperature: 134 °C
- duration of the sterilization hold time: 18 minutes
- Pressure: 2 bars

Storage: Then, store the sterilized products in a dry place protected from dust.

Before re-use, if the integrity of the packaging is not conformable, recondition and re-sterilize in accordance with the defined protocol.

Disposition of the product: Dispose of the product in receptacles for waste materials of healthcare activities involving infectious risks.

2) Before cleaning the handpiece, insert the protective cap in place of the light guide to prevent any liquid entering the handpiece.

3) After each use, check thoroughly that no composite residue has adhered to the light guide. Remove any residue immediately and check that the surface of the light guide has not been altered. If damage is detected, replace the light guide, as the power of the curing light could be reduced significantly.

4) Under normal conditions of use, the power of the curing light does not vary if the battery is properly charged. Consequently there is no need to check the power as in the case of ordinary curing lights.

However, if in doubt, check the power using the light guide.

- 5) The battery unscrews from the handpiece for replacement.
- 6) Remove the light guide to check that the LED is clean and undamaged. Clean it with a dry air jet if necessary.
- 7) Never use an ultrasound bath for cleaning the Mini L.E.D. AutoFocus or its accessories.

## 9 - TROUBLESHOOTING

If a problem occurs, carry out the following preliminary checks:

- 1) Make sure that the base station is correctly connected to the power supply so that the battery recharges normally. If the Mini L.E.D. AutoFocus rest is not lit up even when the connection is correct, contact the after-sales service or your dealer.
- 2) Check the status of the battery charge indicator (4) before pressing the ON/OFF button. If the indicator light is red, an automatic cutout prevents the unit operating. In this case, recharge the battery by placing the curing light on the base station.
- 3) Intensive use of the curing light may cause a rise in temperature inside the unit. If this occurs, an automatic cut-out prevents operation of the unit; the indicator light (4) flashes red, there is an audible beep, and a flashing warning message "OH" is displayed on the LCD screen. Leave the unit to cool down for a few minutes; the curing light should then work normally again.

- 4) If the translucent rest lights up red during the check with the radiometer, clean the light guide and the radiometer window thoroughly to remove any possible traces of material, then try again. Check whether the light guide is intact. If not, it must be replaced.

## 10 - TECHNICAL SERVICE

Your dealer's after-sales service is at your disposal to answer any questions.

## 11 - SPECIFICATIONS

### Unit

Model:	Mini LED AutoFocus
Weight*:	200g
Dimensions*:	Ø24 x 210mm
Classification:	Ordinary
	

\* Without light guide

### Power adapter

Input voltage:	100 - 240 VAC
Frequency:	47 - 63Hz
Output voltage:	12V
Output current:	0.6A
Classification:	II IP 40

### Base station

Supply voltage:	12 VDC
Protection:	Fuse, 2A
Classification:	Continuous operation IPX 0

### Battery

Type:	Lithium-ion
Dimensions:	87 x Ø24mm
Capacity:	2,400 mAh

### Optical specifications

Wavelength:	420-480 nm
Intensity:	2,200mW/cm <sup>2</sup> ±10% (Ø7.5mm guide)
Power:	760mW ±10%
	 Class 2M laser

### Temperatures

Operating:	+10°C to +40°C
Storage:	-20°C to +70°C

### Humidity

Operating:	30% to 75%
Storage:	10% to 100% including condensation

## 12 - ELECTROMAGNETIC COMPATIBILITY

 **Warning:** The charger power cord must be kept apart from those of any nearby devices.

**Mini L.E.D. AutoFocus** requires special precautions to be taken with regard to electromagnetic compatibility. It must be installed and prepared for use as described in chapter 3.

Certain types of mobile telecommunication devices such as mobile telephones are likely to interfere with the **Mini L.E.D. AutoFocus**. The recommended separation distances in this paragraph must therefore be complied with.

The **Mini L.E.D. AutoFocus** must not be used near or on top of another device.

If this cannot be avoided, its operation under the conditions of use must be checked beforehand.

The use of accessories other than those specified or sold by SATELC as replacement parts may have the consequence of increasing the emissions or decreasing the immunity of the **Mini L.E.D. AutoFocus**.

### 12.1 - Electromagnetic emissions

The **Mini L.E.D. AutoFocus** is intended for use in the electromagnetic environment specified in the table below. The user and/or installer must ensure that the **Mini L.E.D. AutoFocus** is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emission - CISPR 11.	Group 1 Class B	<p>Mini L.E.D. AutoFocus uses RF energy for internal operation. Therefore, its radiofrequency emissions are very low and are not likely to cause any interference in nearby equipment.</p> <p>Mini L.E.D. AutoFocus uses RF energy for internal operation. Therefore, its radiofrequency emissions are very low and are not likely to cause any interference in nearby equipment.</p>

## 12.2 - Magnetic and electromagnetic immunity

The Mini L.E.D. AutoFocus is intended for use in the electromagnetic environment specified in the table below.

The user and/or installer must ensure that the Mini L.E.D. AutoFocus is used in such an electromagnetic environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Magnetic field at 50Hz. IEC61000-4-8	3A/m	3A/m	The intensity of the magnetic field should be equivalent to that of a typical commercial or hospital environment (hospital, clinic).
Electrostatic discharge (ESD) IEC 61000-4-2.	± 6kV contact ± 8kV air	± 6kV contact ± 8kV air	Floors must be wood, concrete, cement or tiled. If floors are covered with synthetic material (carpet, etc.), the relative humidity must be at least 30%.
Electrical fast transients IEC 61000-4-4.	± 2kV for power supply lines	± 2kV for power supply lines	Power quality should be that of a typical commercial or hospital environment.
Surges IEC 61000-4-5.	± 1kV differential mode ± 2kV common mode	± 1kV differential mode ± 2kV common mode	Power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations IEC 61000-4-11.	<5% Ur (>95% dip in Ur) for 0.5 cycles. 40% Ur (60% dip in Ur) for 5 cycles 70% Ur (30% dip in Ur) for 25 cycles <5% Ur (>95% dip in Ur) for 250 cycles	<5% Ur (>95% dip in Ur) for 0.5 cycles. 40% Ur (60% dip in Ur) for 5 cycles 70% Ur (30% dip in Ur) for 25 cycles <5% Ur (>95% dip in Ur) for 250 cycles	Power quality should be that of a typical commercial or hospital environment. If the use of the Mini L.E.D. AutoFocus requires continued operation during a power cut, it is recommended that the product be powered from a separate power supply (UPS, etc.).

## 12.3 - Electromagnetic immunity / mobile radiofrequency equipment

The Mini L.E.D. AutoFocus is intended for use in the electromagnetic environment specified in the table below.

The user and/or installer must ensure that the device is used in such an electromagnetic environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Portable and mobile radiofrequency communications devices must not be used near the Mini L.E.D. AutoFocus (including its cables) at a distance less than that recommended and calculated according to the frequency and power of the emitter.			
Conducted disturbance, radiofrequency fields. IEC61000-4-6	3 V/m 150 KHz to 80 MHz	3 V/m	Recommended separation distance: $d = 1.2 \sqrt{P}$
Radiated radiofrequency electromagnetic field. IEC61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80MHz to 800MHz. $d = 2.3 \sqrt{P}$ 800MHz to 2.5GHz. Where P is the maximum power rating of the emitter in watts (W) according to the manufacturer's specifications and d is the recommended minimum separation distance in meters (m).

The electromagnetic field strengths of fixed radiofrequency emitters, as determined by an electromagnetic environment measurement (a), must be less than the compliance level in each frequency range (b).  
Interference may occur near equipment marked with the symbol below:

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These specifications may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and persons.

(a): The electromagnetic field strengths of fixed radiofrequency emitters, such as base stations for mobile telephones (cellular/cordless), mobile radios, amateur radios, AM/FM radio broadcasts and TV broadcasts cannot be determined exactly by theory. To assess the electromagnetic environment due to fixed radiofrequency emitters,

an electromagnetic environment measurement must be made. If the measured radiofrequency field strength in the immediate environment where the product is used exceeds the compliance level specified above, the performance of the product must be tested to verify whether it conforms to the specifications. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the product.

(b): In the 150 kHz to 80 MHz frequency range, the electromagnetic field strengths must be less than 3 V/m.

## 12.4 - Recommended separation distances

The **Mini L.E.D. AutoFocus** is intended for use in an electromagnetic environment in which radiated radiofrequency disturbances are controlled.

The **Mini L.E.D. AutoFocus** user and/or installer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile radiofrequency communications equipment (emitters) and the **Mini L.E.D. AutoFocus**, according to the maximum output power of the equipment, as recommended in the table below.

Rated max. power of the emitter (W)	Separation distance in meters (m) according to emitter frequency		
	150 KHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12 m	0.12 m	0.23 m
0.1	0.38 m	0.38 m	0.73 m
1	1.2 m	1.2 m	2.3 m
10	3.8 m	3.8 m	7.3 m
100	12 m	12 m	23 m

For emitters rated at max. power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the emitter, where P is the max. power rating of the emitter in watts (W) according to the manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These specifications may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and persons.

## 12.5 - Cable lengths

Cables and accessories	Maximum length	Complies with:
charger power cord	<3m	RF emission, CISPR 1 - Class B Immunity to magnetic fields: IEC61000-4-8. Immunity to electrostatic discharge: IEC61000-4-2 Immunity to electrical fast transients/bursts: IEC61000-4-4 Immunity to surges: IEC61000-4-5 Immunity to voltage dips, short interruptions and voltage variations: IEC61000-4-11 Immunity to conducted disturbances induced by radiofrequency fields: IEC61000-4-6 Immunity to radiated radiofrequency electromagnetic fields: IEC61000-4-3

## 13 - DISPOSAL AND RECYCLING

As electrical and electronic equipment, the device must be disposed of according to a specialized procedure for collection, pick-up and recycling or destruction (in particular on the European market, with reference to Directive 2002/96/EC of 23/01/2003).

When your device reaches the end of its life, we consequently recommend that you contact your dental equipment dealer (or, failing this, the nearest ACTEON GROUP office, the list of which is given in chapter 17), for information on how to proceed.

## 14 - LIABILITY

The manufacturer is not liable if:

- the manufacturer's installation recommendations have not been followed (supply voltage, electromagnetic environment, etc.);
- repairs have been performed by persons not authorized by the manufacturer;
- the device has been used in an electrical installation which does not comply with current standards;
- the device has been used in a way which is not stipulated in this Manual;
- accessories other than those supplied by SATELEC have been used;
- the instructions in this document have not been followed.

The manufacturer reserves the right to modify the unit and/or the Operating Manual without notice.

## 15 - REGULATIONS

This medical device is classified as class IIa according to European Directive 93/42/EEC.

This equipment is manufactured in compliance with the current IEC 60601-1 standard.

This equipment has been designed and manufactured according to an ISO 13485-certified quality assurance system.

## 16 - SYMBOLS AND ABBREVIATIONS

SYMBOL	DEFINITION
	Alternating current
	Direct current
	Warning, please refer to the accompanying documentation
	Class 2M laser
	Type B
	Class 2

Note:

Technical personnel of the Satelec authorized dealer network can obtain from **ACTEON Group** on request all the information they need for repair of the parts of the curing light that ACTEON has identified as repairable.

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