m i n i l e d<sup>TM</sup>  
ADVANTAGE

---

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# ENGLISH

# SUMMARY

1- SAFETY INSTRUCTIONS	6
2- PRESENTATION	7
3- INSTALLATION	7
4- PRODUCT DESCRIPTION	7
5- MODES USED	7
6- EVERYDAY USE	8
7- CARE	8
8- TROUBLESHOOTING	8
9- SPECIFICATIONS	9
10- ELECTROMAGNETIC COMPATIBILITY	10
11- DISPOSAL AND RECYCLING	12
12- LIABILITY	12
13- REGULATIONS	12
14- SYMBOLS AND ABBREVIATIONS	13

# 1 - SAFETY INSTRUCTIONS

**Warning:** Before connecting the MINI L.E.D. ADVANTAGE, please read the safety instructions carefully.

The manufacturer's liability for the Mini L.E.D. Advantage is applicable only if it is used as directed and if the safety instructions are properly complied with. Safety instructions mean a good prior knowledge of dentistry and photo-polymerization as well as understanding of the specific guidelines given in this User's Manual.

- 1) This device must be used in accordance with the instructions contained in this user's manual. We do not accept liability in the case of non-compliance with the instructions or if the device is used for any purpose other than that recommended.
- 2) Before plugging in the device the base station to the mains, check that the mains voltage is set between the minimum and the maximum as indicated on the label on the mains adapter. A voltage outside the range stated will damage the device and could injure the operator and/or the patient.
- 3) The light rays produced by the device can be harmful and should never be directed at the eyes even if adequate protective glasses are worn (Laser class 2M). The light is only intended for treatment inside the mouth.
- 4) Condensation inside an electrical appliance can be dangerous. If the device needs to be transported from a cold to a warm area, it should not be used immediately but should be left long enough to reach the new room temperature.
- 5) Unless otherwise stated in this manual, do not introduce or try to introduce any metal objects into the device or the base station. This could cause an electric shock.
- 6) The manufacturer cannot be held liable if damaged parts or accessories are replaced by parts other than those specified by the manufacturer. In particular, use of any other light guide, mains adaptor or batteries other than those originally supplied could be dangerous for the patient/user and/or for the correct use of the instrument.
- 7) In the case of a fault, unplug the device immediately and ensure that no other person can use it before it has been checked by the manufacturer or the manufacturer's authorized dealer. A fault could be due to incorrect use or physical damage.
- 8) Do not use the MINI L.E.D. ADVANTAGE close to any heat source. Use of solvent or flammable products on this device could cause damage.
- 9) The use of detergents on the device for cleaning or any other purposes may cause a short circuit.
- 10) If repairs are necessary, the device should only be opened by the manufacturer's authorized dealer.
- 11) Any device used for polymerization, including the MINI L.E.D. ORTHO, must not be used on, or by, any person who suffers, or has suffered in the past, from photo-biological reactions (including urticaria solaris or protoporphyrine erythropoietique), or any person receiving medical treatment and using photosensitive medication (including methoxsalenes or dymethyl-chlorotetracycline).
- 12) Any person, practitioner or patient, who has previously suffered from a retina or crystalline condition or who has undergone eye surgery, in particular cataract treatment, should visit an eye specialist before using the device or before using the MINI L.E.D. ADVANTAGE. Even with the patient's agreement, great care should be taken because the intensity of the light is potentially dangerous. It is strongly recommended that protective glasses suited to using equipment emitting class 2M type laser rays with a wave length of between 420 and 480nm are used at all times.
- 13) Any variation in the voltage of the power supply or in the electromagnetic field outside the normal safety regulations could cause the device to come on automatically, or interfere with its operation.
- 14) Electromagnetic interference: the device has been tested and approved for all the current standards (EN 60 601-1-2) both with regard to emission (electrical field rays and supply voltage interference) and immunity (electrical fields rays, intermittent bursts on the supply, electrostatic discharges and supply shock waves).
- 15) MINI L.E.D. ADVANTAGE must not be used if the patient and/or the operator has a cardiac stimulator or any other active implant (e.g. a cochlear implant).
- 16) The device is not designed to withstand shocks

delivered by an electric defibrillator.

- 17) If required, your dealer can supply all spare parts.  
Please note that the battery should always be returned to the dealer.
- 18) For use only by a dental professional.
- 19) Make sure that the battery is unscrewed before transporting the device.

## 2 - PRESENTATION

The **MINI L.E.D. ADVANTAGE** produces a visible blue light between a wavelength of 420 and 480nm for the photo-polymerization of dental materials. After checking the packaging for damage (please retain the packaging), you will find the following:

- the handpiece **MINI L.E.D. ADVANTAGE**, fitted with a protective cap, incorporating the light-emitting diode,
- The power supply cable, consisting of the mains adapter and the power socket.
- The multifibre light guide, opalescent, curvature 45°, Ø7.5mm, sterilizable
- The protective light shield
- The stand for the handpiece
- The User's Manual and other documents.

### Available options:

- a multifibre light guide, opalescent, curvature 45°, Ø5.5mm, sterilizable

## 3 - INSTALLATION

- 1- Prior to each use, sterilize the light guide and disinfect both base station and handpiece (see also the chapter on care).
- 2- Screw the power supply cable into the handpiece and insert the sterilized light guide. Ensure that the light guide is inserted correctly (the friction contact makes a click).
- 3- Plug in the power supply cable to the mains supply, after checking its compatibility with your surgery's power supply voltage.

The handpiece's green LED lights up and the device beeps 3 times.

## 4 - PRODUCT DESCRIPTION

Now that the **MINI L.E.D. ADVANTAGE** is fully charged, you need to familiarize yourself with its buttons and menus:

### Handpiece

#### a) LED indicators

The handpiece has four LED indicators:

- One LED located at the top, next to the light guide, indicates the status:

Color	Mode
Green	Normal mode
Red	Power supply problem
Red flashing	Thermal protection

- A row of three LEDs indicates the curing mode selected

Symbol	Mode
	Fast
	Pulse
	Stepped

#### b) Buttons

The handpiece has two buttons:

- one navigation button which is used to switch from one menu to another
- one on/off light button (trigger) close to the light guide.

## 5 - MODES USED

The **MINI L.E.D. ADVANTAGE** has three different modes. To select the mode you want to use, press the navigation button. The green LED indicator next to each symbol shows the mode selected.

**"FAST CURING" mode:** when the handpiece is removed from the base station, the device will be set on fast curing program. This menu activates maximum light power (depending on the light guide) for 10 seconds.

**"PULSE CURING" mode:** Selected by means of the navigation key. This menu gives full power in a pulsed

mode, which consists of an emission of 10 successive one second flashes with a rest period of 250ms between the flashes.

**"STEPPED CURING" mode:** Selected by means of the navigation key. "Stepped curing" (curing similar to the step menu in a halogen device), gives you:

- a soft start of 10 seconds
- full power for 10 seconds

## 6 - EVERYDAY USE

- 1) Each time it is plugged in, the **MINI L.E.D. ADVANTAGE** will be activated in the fast curing mode. For further details, please refer to chapter 5.
- 2) Your **MINI L.E.D. ADVANTAGE** is ready to use once the menu has been selected.
- 3) Place the light guide as close as possible to the surface of the VLC material without touching it, as this could greatly damage the light guide and decrease its efficiency (staining)
- 4) Press the On/Off button briefly to start the curing cycle) The setting will be confirmed by an audible signal (beep).
- 5) At the end of the cycle, you will hear another beep. However, the cycle can be interrupted at any time simply by pressing the On/Off button.
- 6) If the device is not used, after three minutes it goes into Standby mode and all the indicator lights go out. Press any button to reactivate the device.

## 7 - CARE

**Warning:** Disconnect the **MINI L.E.D. ADVANTAGE** properly before using a disinfectant. Clean and sterilize the light guide between patients.

- 1) The **MINI L.E.D. ADVANTAGE** and its accessories can be cleaned and disinfected with ready-to-use cleaning and disinfecting wipes impregnated with alcohol, an amphoteric and a biguanide (refer to manufacturer's instructions) for at least two minutes. Leave the product to take effect for at least 15 minutes. Use wipes that have CE marking or that are in compliance with any standard required by the national regulations.

Wipe dry using a disposable clean non-woven cloth until there are no traces of liquid.

Packaging: Keep sterilizable accessories in single-use sterilization sleeves or bags in compliance with the specifications defined in standard EN ISO 11607-1 or in any equivalent standard required by the national regulations.

Sterilize the light guide and the flexible optical shield individually in an autoclave before each patient, according to the following parameters:

- autoclave, compliant with standard EN 13060, class B
- sterilization temperature: 134 °C
- duration of the sterilization hold time: 18 minutes
- Pressure: 2 bars

Storage: Then, store the sterilized products in a dry place protected from dust.

Before re-use, if the integrity of the packaging is not conformable, recondition and re-sterilize in accordance with the defined protocol.

Disposition of the product: Dispose of the product in receptacles for waste materials of healthcare activities involving infectious risks.

**2) Before cleaning the handpiece, insert the protective cap in place of the light guide to prevent any liquid entering the handpiece.**

3) After each use, check thoroughly that no composite residue has adhered to the light guide. Remove any residue immediately and check that the surface of the light guide has not been altered. If damage is detected, replace the light guide, as the power of the curing light could be reduced significantly.

4) Remove the light guide and check that the LED is clean and undamaged. Clean it with a dry air jet if necessary.

5) Never apply ultrasonic cleaning on the **MINI L.E.D. ADVANTAGE** and its accessories.

## 8 - TROUBLESHOOTING

If a problem appears, before calling your dealer's after sales service:

- 1) Check that the connecting cable is properly connected to the mains power supply. If the LED is not lit up in spite of proper connection, contact the after sales service.
- 2) Check that the status LED shows that the device is properly powered (green light) before pressing the

On/Off button. If the LED is red, an automatic safety mechanism will stop the device from being activated. In this case, check the module's mains power connection and also that the connecting cable is properly screwed into the handpiece.

- 3) Intensive use of the lamp will generally cause a rise in its external temperature that can be felt by the hand holding it. This is absolutely normal. If the temperature inside the **MINI L.E.D. ADVANTAGE** should rise too high, an automatic safety mechanism will stop the lamp, and the status indicator will flash red. Leave the device for a few minutes to cool down and then carry on as normal.
- 4) If the curing is insufficient, check that the LED and the light guide are clean. Use a jet of dry air to remove any dust. If the problem persists or if the light guide is damaged, the **MINI L.E.D. ADVANTAGE** must be returned to the after-sales service department.

Your dealer's after sales service is at your disposal to help with any technical problem.

## 9 - SPECIFICATIONS

### Handpiece

Model:	<b>MINI L.E.D. ADVANTAGE</b>
Weight:	160 g
Size:	Ø 23 x 180 mm
Classification:	IPX0

### Power supply cable

Brand:	Satelec
Model:	Advantage Power supply (Ref. F02459)

### Input:

Voltage:	100 - 240 VAC ~
Current:	400 mA
Frequency:	50-60 Hz

### Output:

Voltage:	5 V DC
Current:	2,4 A
Classification :	Ordinary Class II IP41 Direct current

### Optical specifications

Wavelength:	420-480nm
Intensity:	1250 to 2000 mW/cm <sup>2</sup> +10%
(depending on light guide)	

Class 2M laser

### Temperatures:

Operating:	+10°C to +40°C
Storage:	-20°C to +70°C

### Humidity:

Operating:	30% to 75%
Storage:	10% to 100% including condensation

## 10 - ELECTROMAGNETIC COMPATIBILITY

**Important:** **MINI L.E.D. ADVANTAGE** conforms to the electromagnetic compatibility standards currently in force within the following configuration:

- Ferrite magnet - see wiring instructions.
- Conversion module
- Length of **MINI L.E.D. ADVANTAGE** cord = max. 3m.
- Power supply voltage : 24VAC.
- **MINI L.E.D. ADVANTAGE**

**MINI L.E.D. ADVANTAGE** requires special precautions to be taken with regard to electromagnetic compatibility. It must be installed and prepared for use as described in chapter 3.

Certain types of mobile telecommunication devices such as mobile telephones are likely to interfere with the **MINI L.E.D. ADVANTAGE**.

The recommended separation distances in this paragraph must therefore be complied with.

The **MINI L.E.D. ADVANTAGE** must not be used near or on top of another device. If this cannot be avoided, its operation under the conditions of use must be checked beforehand.

The use of accessories other than those specified or sold by SATELEC as replacement parts may have the consequence of increasing the emissions or decreasing the immunity of the **MINI L.E.D. ADVANTAGE**.

### 10.1 - Electromagnetic emissions

The **MINI L.E.D. ADVANTAGE** is intended for use in the electromagnetic environment specified in the table below. The user and/or installer must ensure that the **MINI L.E.D. ADVANTAGE** is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emission - CISPR 11.	Group 1 Class B	<p><b>MINI L.E.D. ADVANTAGE</b> uses RF energy for internal operation. Therefore, its radiofrequency emissions are very low and are not likely to cause any interference in nearby equipment.</p> <p><b>MINI L.E.D. ADVANTAGE</b> is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.</p>

### 10.2 - Magnetic and Electromagnetic immunity

The **MINI L.E.D. ADVANTAGE** is intended for use in the electromagnetic environment specified in the table below. The user and/or installer must ensure that the device is used in such an electromagnetic environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Magnetic field at 50Hz. IEC61000-4-8	3A/m	3A/m	The intensity of the magnetic field should be equivalent to that of a typical commercial or hospital environment (hospital, clinic).
Electrostatic discharge (ESD) IEC 61000-4-2.	± 6kV contact ± 8kV air	± 6kV contact ± 8kV air	Floors must be wood, concrete, cement or tiled. If floors are covered with synthetic material (carpet, etc.), the relative humidity must be at least 30%.
Electrical fast transients IEC 61000-4-4.	± 2kV for power supply lines	± 2kV for power supply lines	Power quality should be that of a typical commercial or hospital environment.
Surges IEC 61000-4-5.	± 1kV differential mode ± 2kV common mode	± 1kV differential mode ± 2kV common mode	Power quality should be that of a typical commercial or hospital environment.

## 10.3 - Electromagnetic immunity / mobile radiofrequency equipment

**MINI L.E.D. ADVANTAGE** is intended for use in the electromagnetic environment specified in the table below. The user and/or installer must ensure that the **MINI L.E.D. ADVANTAGE** is used in such an electromagnetic environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Portable and mobile radiofrequency communications devices must not be used near the <b>MINI L.E.D. ADVANTAGE</b> (including its cables) at a distance less than that recommended and calculated according to the frequency and power of the emitter.			
Conducted disturbance, radiofrequency fields. IEC61000-4-6	3 V/m 150 KHz to 80 MHz	3 V/m	Recommended separation distance: $d = 1.2 \sqrt{P}$
Radiated radiofrequency electromagnetic field. IEC61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80MHz to 800MHz. $d = 2.3 \sqrt{P}$ 800MHz to 2.5GHz. Where P is the maximum power rating of the emitter in watts (W) according to the manufacturer's specifications and d is the recommended minimum separation distance in meters (m).
The electromagnetic field strengths of fixed radiofrequency emitters, as determined by an electromagnetic environment measurement (a), must be less than the compliance level in each frequency range (b). Interference may occur near equipment marked with the symbol below:			
			

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These specifications may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and persons.

(a): The electromagnetic field strengths of fixed radiofrequency emitters, such as base stations for mobile telephones (cellular/cordless), mobile radios, amateur radios, AM/FM radio broadcasts and TV broadcasts cannot be determined exactly by theory. To assess the electromagnetic environment due to fixed radiofrequency emitters, an electromagnetic environment measurement must be made. If the measured radiofrequency field strength in the immediate environment where the product is used exceeds the compliance level specified above, the performance of the product must be tested to verify whether it conforms to the specifications. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the product.

(b): In the 150 kHz to 80 MHz frequency range, the electromagnetic field strengths must be less than 3 V/m.

## 10.4 - Recommended separation distances

The **MINI L.E.D. ADVANTAGE** is intended for use in an electromagnetic environment in which radiated radiofrequency disturbances are controlled.

**MINI L.E.D. ADVANTAGE** user and/or installer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile radiofrequency communications equipment (emitters) and the **MINI L.E.D. ADVANTAGE**, according to the maximum output power of the equipment, as recommended in the table below.

Rated max. power of the emitter (W)	Separation distance in meters (m) according to emitter frequency		
	To 150 KHz to 80 MHz $d = 1.2 \sqrt{P}$	To 80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	To 800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12 m	0.12 m	0.23 m
0.1	0.38 m	0.38 m	0.73 m
1	1.2 m	1.2 m	2.3 m
10	3.8 m	3.8 m	7.3 m
100	12 m	12 m	23 m

For emitters rated at max. power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the emitter, where P is the max. power rating of the emitter in watts (W) according to the manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These specifications may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and persons.

## 10. 5 Cable lengths

Cables and accessories	Maximum length	Complies with:
Handpiece cord	<3m	RF emission, CISPR 1 - Class B Immunity to magnetic fields: IEC61000-4-8. Immunity to electrostatic discharge: IEC61000-4-2 Immunity to electrical fast transients/bursts: IEC61000-4-4 Immunity to surges: IEC61000-4-5 Immunity to voltage dips, short interruptions and voltage variations: IEC61000-4-11 Immunity to conducted disturbances induced by radiofrequency fields: IEC61000-4-6 Immunity to radiated radiofrequency electromagnetic fields: IEC61000-4-3

## 11 - DISPOSAL AND RECYCLING

As electrical and electronic equipment, the device must be disposed of according to a specialized procedure for collection, pick-up and recycling or destruction (in particular on the European market, with reference to Directive 2002/96/EC of 23/01/2003).

When your device reaches the end of its life, we consequently recommend that you contact your dental equipment dealer (or, failing this, the nearest ACTEON GROUP office, the list of which is given in chapter 15), for information on how to proceed.

## 12 - LIABILITY

The manufacturer is not liable if:

- the manufacturer's installation recommendations have not been followed (supply voltage, electromagnetic environment, etc.);
- repairs have been performed by persons not authorized by the manufacturer;
- the device has been used in an electrical installation which does not comply with current standards;
- the device has been used in a way which is not stipulated in this Manual;
- accessories other than those supplied by SATELEC have been used;
- the instructions in this document have not been followed.

The manufacturer reserves the right to modify the unit and/or the Operating Manual without notice.

## 13 - REGULATIONS

This medical device is classified as class IIa according to European Directive 93/42/EEC.

This equipment is manufactured in compliance with the current IEC 60601-1 standard.

This equipment has been designed and manufactured according to an ISO 13485-certified quality assurance system.

## 14 - SYMBOLS AND ABBREVIATIONS

SYMBOL	DEFINITION
	Alternating current
	Direct current
	Warning, please refer to the accompanying documentation
	Device emitting Class 2M laser rays
	Type B
	Class II

### Note:

Technical personnel of the **ACTEON Group** authorized dealer network can obtain from SATELEC on request all the information they need for repair of the parts of the equipment that ACTEON has identified as repairable.

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